Алматы (7273)495-231 Ангарск (3955)60-70-56 Архангельск (8182)63-90-72 Астрахань (8512)99-46-04 Барнаул (3852)73-04-60 Белгород (4722)40-23-64 Благовещенск (4162)22-76-07 Брянск (4832)59-03-52 Владивосток (423)249-28-31 Владикавказ (8672)28-90-48 Владикавказ (8672)28-90-48 Волоград (844)278-03-48 Вологда (8172)26-41-59 Воронеж (473)204-51-73 Екатеринбург (343)384-55-89 Ижевск (3412)26-03-58 Иваново (4932)77-34-06 Иркутск (395)279-98-46 Казань (843)206-01-48 Калининград (4012)72-03-81 Капуга (4842)92-23-67 Кемерово (3842)65-04-62 Киров (8332)68-02-04 Коломна (4966)23-41-49 Кострома (4942)77-07-48 Краснодар (661)203-40-90 Красноярск (391)204-63-61 Курск (4712)77-13-04 Курган (352)50-90-47 Липецк (4742)52-20-81

Киргизия (996)312-96-26-47

Магнитогорск (3519)55-03-13 Москва (495)268-04-70 Мурманск (8152)59-64-93 Набережные Челны (8552)20-53-41 Нижний Новгород (831)429-08-12 Новокузнецк (3843)20-46-81 Ноябрьск (3496)41-32-12 Новосибирск (383)227-86-73 Ноябрьск (3496)41-32-12 Омск (3812)21-46-40 Орел (4862)44-53-42 Оренбург (3532)37-68-04 Пенза (8412)22-31-16 Петрозаводск (8142)55-98-37 Псков (8112)59-10-37

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Сургут (3462)77-98-35
Тамбов (4752)50-40-97

Тверь (4822)63-31-35 Тольяти (8482)98-41-53 Тула (4872)98-41-53 Тула (4872)33-79-87 Тюмень (3452)66-21-18 Улан-Уда (3012)59-97-51 Ульяновск (8422)24-23-59 Уфа (347)229-48-12 Хабаровск (4212)92-98-04 Чебоксары (8352)28-53-07 Черяповец (8202)49-0-64 Чита (3022)38-34-83 Якутск (4112)23-90-97 Ярославль (4852)69-52-93

https://schmitz.nt-rt.ru || szk@nt-rt.ru

# Тележки Varimed ISO

Россия (495)268-04-70



# varimed® ISO trolleys Our experts for ward logistics.

SCHMITZ's market-leading ISO trolleys offer the perfect assistance for everyday clinical practice. Hardwearing and flexible – varimed® ISO trolleys are used in almost every German clinic.

# **Optimal organisation**

A customised range of options ensures systems perfectly tailored to your structures – we organise your workplace.

# High flexibility

Different heights, different drawer configurations, full-size or half-size ISO modules – the SCHMITZ ISO trolley adapts to your daily work processes.

# Hygiene taken a step further

Removable drawer panels and closed surfaces enable rapid, complete and simple cleaning – a cleaning standard redefined by SCHMITZ.

# Perfect hygiene standard

Nearly complete elimination of germs thanks to antimicrobial coatings – no chance for multidrug-resistant pathogens.

# As individual as your requirements

Our colour and configuration concept guarantees ideally matched design options.



















# varimed® ISO trolleys

# Wall protecting wheels protect the trolleys, furniture and the walls from damage

Closed double easy-roll ball-bearing castors, non-contact marking and abrasion-free quality. Easy-to-manoeuvre, silent ISO trolleys.

## Side rails (202.0016.0)

for easy attachment of accessories.

# Hidden drawer slides and full pullout function for optimal hygiene

Drawers with automatic pull-in function and Airmatic shock absorber system. Easy to remove front makes cleaning simple

# **Ergonomic guiding bar in two models:**

- Straight guiding bar (as shown on the left)
- Curved model

The guiding bar can be fitted to the right or left side or can be left off.

### Germs stand no chance


Our special antimicrobial coating fights any local germs effectively, leaving them no place to spread. For more information on antimicrobial coating, see pages 168 and 169.

Surrounding silicone sealing profile, disinfectant-proof Opening angle of the doors max. 125°.

### ISO rail walls

allow for flexible placement of SCHMITZ ISO modules.

















# varimed® ISO trolleys

## Practical design

The drawers can be fitted with one full-size ISO module or two half-size ISO modules.

# ISO trolley models with a height of 850 mm or 1,070 mm, tailored to your processes

For unequalled flexibility: ISO trolleys 850 mm high can be fitted with up to 8 full-size or 16 half-size ISO modules; ISO trolleys 1,070 mm high can be fitted with up to 11 full-size or 22 half-size ISO modules.

An overview of the drawer layout options for full-size and half-size ISO modules can be found on page 167.

## Easy exchange of modules

for fast replenishment of ISO trolleys. The modules can be removed without tilting or tipping.

All ISO module drawers have recessed grips, making it much easier to remove the modules.

# Take safety and security even further!

Innovative lock designs with state-of-the-art technology.



# varimed® ISO trolleys Width 720 mm, height 850 mm

Model No. Description

### 260.0025.0 ISO Trolley

with 4 x drawers for modules of max. 100 mm height. 4 closed double castors Ø 125 mm, plastic, ball bearing, 2 of which electroconductive, 2 full lock-type castors, and 1 with directional locking, silver grey frame cover with revolving bumpers, plastic top (silver grey) with profile edge, guiding bar.

Modules not included.



### 260.0026.0 ISO Trolley

as ISO trolley 260.0025.0, however with 2 x drawers for modules of max. 100 mm height and 1 x drawer for modules of max. 200 mm height.

Modules not included.



### 260.0027.0 ISO Trolley

as ISO trolley 260.0025.0, however with 2 x drawers for modules of max. 50 mm height and 3 x drawers for modules of max. 100 mm height.









# varimed® **ISO** trolleys Width 720 mm, height 850 mm

Model No. Description

#### 260.0028.0 **ISO Trolley**

as ISO trolley 260.0025.0, however with 2 x drawers for modules of max. 50 mm height, 1 x drawer for modules of max. 100 mm height and 1 x drawer for modules of max. 200 mm height.

Modules not included.

#### 260.0029.0 **ISO Trolley**

as ISO trolley 260.0025.0, however with 2 x drawers for modules of max. 100 mm height and 4 x drawers for half modules of max. 100 mm height.

Modules not included.

#### 260.0030.0 **ISO Trolley**

as ISO trolley 260.0025.0, however with 2 x drawers for modules of max. 100 mm height and 2 x drawers for half modules of max. 200 mm height.

Modules not included.

#### **ISO Trolley** 260.0038.0

as ISO trolley 260.0025.0, however with 2 hinged doors. Prepared for installation of ABS rail walls for full-size ISO modules (e.g. model no. 202.0082.0 from SCHMITZ or those of another manufacturer). We recommend the use of stops 202.0072.0 (pair) for full-size ISO modules as accessory.









# varimed® ISO trolleys Width 720 mm, height 1,070 mm

Model No. Description

### 260.0031.0 ISO Trolley

with 1 x drawer for modules of max. 50 mm height and 5 x drawers for modules of max. 100 mm height. 4 closed double castors Ø 125 mm, plastic, ball bearing, 2 of which electroconductive, 2 full lock-type castors, and 1 with directional locking, silver grey frame cover with revolving bumpers, plastic top (silver grey) with profile edge, guiding bar. **Modules not included.** 



### 260.0032.0 ISO Trolley

as ISO trolley 260.0031.0, however with 1 x drawer for modules of max. 50 mm height, 3 x drawers for modules of max. 100 mm height and 1 x drawer for modules of max. 200 mm height.

Modules not included.



## 260.0033.0 ISO Trolley

as ISO trolley 260.0031.0, however with 3 x drawers for modules of max. 50 mm height and 4 x drawers for modules of max. 100 mm height.



# varimed® ISO trolleys Width 720 mm, height 1,070 mm

Model No. Description

# 260.0034.0 ISO Trolley

as ISO trolley 260.0031.0, however with 3 x drawers for modules of max. 50 mm height, 2 x drawers for modules of max. 100 mm height and 1 x drawer for modules of max. 200 mm height.


Modules not included.



### 260.0035.0 ISO Trolley

as ISO trolley 260.0031.0, however with 1 x drawer for modules of max. 50 mm height, 1 x drawer for modules of max. 100 mm height and 8 x drawers for half modules of max. 100 mm height.





## Model No. Description

# 260.0036.0 ISO Trolley

as ISO trolley 260.0031.0, however with 1 x drawer for modules of max. 50 mm height, 2 x drawers for modules of max. 100 mm height, 2 x drawers for half modules of max. 100 mm height and 2 x drawers for half modules of max. 200 mm height.

Modules not included.

### 260.0037.0 ISO Trolley

as ISO trolley 260.0031.0, however with 2 hinged doors. Prepared for installation of ABS rail walls for full-size ISO modules (e.g. model no. 202.0086.0 from SCHMITZ or those of another manufacturer). We recommend the use of stops 202.0072.0 (pair) for full-size ISO modules as accessory.

Modules not included.

# varimed® ISO trolleys Width 720 mm, height 1,070 mm

Model No. Description

# 260.0050.0 ISO Trolley

as ISO trolley 260.0031.0, however with 5 x drawers for modules of max. 100 mm height, lower drawer with raised front panel.

Modules not included.



# 260.0051.0 ISO Trolley

as ISO trolley 260.0031.0, 3 x drawers for modules of max. 100 mm height, of which the 3<sup>rd</sup> with raised front panel, and 1 x drawer for modules of max. 200 mm height.

Modules not included.



# varimed® ISO trolleys Width 360 mm, height 1,070 mm

Model No. Description

### 260.0052.0 ISO Trolley

as ISO trolley 260.0031.0, however, **360 mm** body width, with 5 x drawers for half modules of max. 100 mm height, lower drawer with raised front panel.

Modules not included.



## 260.0039.0 ISO Trolley

as ISO trolley 260.0031.0, however, **360 mm body width**, with 1x drawer for half modules of max. 50 mm height, 3x drawers for half modules of max. 100 mm height and 1x drawer for half modules of max. 200 mm height.

Modules not included.



Алматы (7273)495-231 Ангарск (3955)60-70-56 Архангельск (8182)63-90-72 Астрахань (8512)99-46-04 Барнаул (3852)73-04-60 Белгород (4722)40-23-64 Благовещенск (4162)22-76-07 Брянск (4832)59-03-52 Владивосток (423)249-28-31 Владикавказ (8672)28-90-48 Владикавказ (8672)28-90-48 Вологорад (844)278-03-48 Вологорад (8172)26-41-59 Воронеж (473)204-51-73 Екатеринбург (343)384-55-89 Ижевск (3412)26-03-58 Иваново (4932)77-34-06 Иркутск (395)279-98-46 Казань (843)206-01-48 Калининград (4012)72-03-81 Кануга (4842)92-23-67 Кемерово (3842)65-04-62 Киров (8332)68-02-04 Коломна (4966)23-41-49 Кострома (4942)77-07-48 Краснодар (861)203-40-90 Курасноярск (391)204-63-61 Курск (4712)77-13-04 Курган (3522)50-90-47 Липецк (4742)52-20-81

Киргизия (996)312-96-26-47

Москва (495):268-04-70 Мурманск (8152):59-64-93 Набережные Челны (8552):20-53-41 Нижний Новгород (831):429-08-12 Новокузнецк (3843):20-46-81 Ноябрьск (3496):41-32-12 Новосибирск (383):227-86-73 Ноябрьск (3496):41-32-12 Омск (3812):21-46-40 Орел (4862):44-53-42 Оренбург (3532):37-68-04 Пенза (8412):22-31-16 Петрозаводск (8142):55-98-37 Псков (8112):59-10-37

Магнитогорск (3519)55-03-13

Россия (495)268-04-70

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Сургут (3462)77-98-35
Тамбов (4752)50-40-97

**К**азахстан (772)734-952-31

Тверь (4822)63-31-35 Тольяти (8482)63-91-07 Томск (3822)98-41-53 Тула (4872)33-79-87 Тюмень (3452)66-21-18 Улан-Удэ (3012)59-97-51 Ульяновск (8422)24-23-59 Уфа (347)229-48-12 Хабаровск (4212)92-98-04 Чебоксары (8352)28-53-07 Челябинск (351)20-20-61 Череповец (8202)49-02-64 Чита (3022)38-34-83 Якутск (4112)23-90-97 Ярославль (4852)69-52-93